

Kanalas, csöves széklet-mintavevő használata

A mintavétel után 24 órán belül kérjük a mintát a laboratóriumban leadni.

A minta addig szobahőn tárolható. Hasmenéses széklet nem alkalmas a vizsgálatra. Menstruáció alatt ne történjen mintavétel!



SZÉKLET PANCREAS SPECIFIKUS ELASZTÁZ ÉS ZONULIN MINTAVÉTEL

A mintavétel menete:

1. Csavarja le a cső a kupakját!
2. A műanyag pálca kanalas végével kanalazzon annyi székletmintát a csőbe, hogy az a kanál feletti műanyag peremig érjen.
3. Zárja vissza a csövet!
4. Mintavétel után a lehető legrövidebb időn belül juttassa vissza a mintát a laboratóriumba! Ha ez nem lehetséges, akkor a mintát fagyasztani szükséges.

Use of the stool sampling spatula

The sample should be returned to one of our clinics within 24 hours.

The sample can be stored at room temperature for 24 hours. Diarrheal stool is inadequate for the test. Sample collection should not be done during menstrual period.



STOOL SAMPLING FOR PANCREATIC ELASTASE AND ZONULIN

Sample collection:

1. Open the tube by removing the cap.
2. Use the plastic spatula to place stool inside the glass tube. Its level should reach the rim over the spoon (walnut size).
3. Close the tube by tightening the cap.
4. After sample collection, please return the sample to one of our clinics as soon as possible. In case the sample cannot be returned within 24 hours after the sample collection, please store the sample frozen (-20°C).